

Typenex® TypeSafe Segment Sampling Device

transfusion process steps:

Patient
Identification

Sample
Identification

Blood Bank
Component

Blood
Unit Issue

Recipient
Verification

Lab Supply
Basics

- TypeSafe™ allows for segment puncture without scissors for cross-matching procedures in the lab
- The large cannula allows cells to pass through with ease
- Top and bottom of device are easily identified by technicians in the blood bank
- Ribbed exterior allows for easy handling with gloved hands
- Fits into 10 and 12mL tube sizes
- Wide mouth fits most standard segment sizes

wide
mouth


ribbed
exterior

tapered
mid-section

large
cannula



product is available in these colors:



**4R4750
BLUE**

- **Highlights: FDA Cleared**
- **Indications: Single, In-Vitro Diagnostic Use**
- **Packaging: 2 bags/500 units each,
1,000 units/case**

Enable your unique transfusion process with a Typenex blood band accessory solution. Contact us for more information.

Phone | 866 897 3639

Fax | 866 897 3630

www.typenex.com

AN-0703-014.1